***njphillips***<sup>™</sup>

**DATAMARS**

**65mL FIXED DOSE POUR-ON WITH BENT  
NOZZLE (2.5mL INCREMENTS)**

**INSTRUCTION AND TROUBLESHOOTING GUIDE**

**ENGLISH | FRANÇAIS | ITALIANO | DEUTSCH | SVENSKA | PORTUGUÊS**

# 65ML FIXED DOSE, 2.5ML INCREMENTS POUR-ON WITH BENT NOZZLE


**This instrument is designed for external application only.**

As components in this instrument may be affected by solvents in other formulations, no responsibility will be accepted by the manufacturer should the instrument be used with products other than those intended.

## INSTRUCTIONS FOR USE:

### 1 CONNECTIONS

Attach the enclosed connecting tube to the 30ml applicator and to the draw off system of the container. Make sure the springs provided are screwed over the tubing in both locations in a counter-clockwise direction. This will prevent the tubing from kinking at these connection points.

### 2 PRIMING THE APPLICATOR


Once the applicator is connected to the container, it **MUST** be primed on maximum dose to ensure all air is expelled from the cylinder. Set the dose adjustor ring to 30ml and depress the lever until the formulation completely fills the applicator cylinder. **The applicator must be held vertically, with the nozzle pointed upwards, to ensure the applicator is fully primed.**

Care must be taken to ensure the liquid does not come into contact with any part of the operators body. Chemicals may cause injury to the operator. Once the applicator is primed, then change the dose as required. **DO NOT** squeeze the lever to change the dose. For more information see MSDS.

### 3 APPLICATION

Refer to the chemical manufacturer's label and instructions for the correct application position and dose recommendations. The illustration shown, indicates the preferred position of a 5

Litre upright draw-off system and 20 Litre container.


#### LEGEND:

- |                               |                              |                            |
|-------------------------------|------------------------------|----------------------------|
| 1. Fan Nozzle                 | 11. Piston                   | Ring                       |
| 2. Diffused Nozzle Tip        | 12. Piston Washer            | 23. Inlet Adaptor          |
| 3. Delivery Valve Spring      | 13. Dose Adjustor            | 24. Lever Pads             |
| 4. Delivery Valve & Seal Ring | 14. Dose Sleeve              | 25. Return Spring          |
| 5. Delivery Cage Lock Nut     | 15. Dose Adjustor Guide Ring | 26. Lever Pin              |
| 6. Nozzle Tip Seal Ring       | 16. Push Rod Seal Ring       | 27. Lever                  |
| 7. Cylinder                   | 17. Push Rod                 | 28. Barb Tube Adaptor      |
| 8. Cylinder Seal Ring         | 18. Inlet Adaptor Seal Ring  | 29. Container Tube Adaptor |
| 9. Nozzle                     | 19. Handle                   | 30. Feed Tube Spring       |
| 10. Piston Washer             | 20. Dose Collar              |                            |
|                               | 21. Inlet Valve Spring       |                            |
|                               | 22. Inlet Valve & Seal       |                            |

**1** As the dose is measured from the rear of the cylinder, some of the formulation will be visible on completion of delivery stroke for doses set at less than the maximum 65ml dose.

**2** The cylinder of this instrument is removed by unscrewing. Should the cylinder be removed, care **must** be taken when replacing the cylinder to ensure the piston washers are not folded back. Incorrect assembly of the piston washers into the cylinder will affect the applicators functionality. See diagram below.


### CARE AND MAINTENANCE:




Attention to cleanliness is essential to ensure continued high performance from this instrument. After each use, or after storage for long periods between use, maintain the applicator by flushing thoroughly with a soap and water mix, then water only. If at any time the instrument should become sluggish in operation, maintenance by cleaning and lubrication with NJ Phillips Lubricant should overcome the problem.

**DO NOT** store your applicator or feedtube full of product. Clean as per the "Care and Maintenance" instructions.

## QUICK REFERENCE TROUBLESHOOTING GUIDE

SYMPTOM	PROBABLE CAUSE	CORRECTIVE ACTION
1a. Product being returned to the container from the applicator.	Foreign matter lodged under the inlet valve.	Remove cylinder (7) and inspect inlet valve to ensure it is clean. Lubricate with NJ Phillips Lubricant and reassemble cylinder.
	Foreign matter lodged under delivery valve.	Remove the nozzle (1 or 2) and ensure the delivery holes are clean and clear. Clean the delivery valve (4) and spring (3) and ensure it is positioned in the correct direction.
1b. Unable to draw product from the container.	Kink in feed tube or blockage within the tube or draw-off system.	Ensure the feed tube spring is located correctly on the applicator and draw-off. 
2. Product leaking out of the nozzle or air being drawn into the cylinder from the nozzle end.	Nozzle seal ring is damaged.	Replace the nozzle seal ring (6).
	Delivery valve and spring are incorrectly assembled.	Reassemble the valve (4) and spring (3) correctly (see parts illustration).
	Delivery valve sealing edge damaged.	Replace the delivery valve (4) and spring (3).
	Delivery cage seal ring damaged.	Replace the delivery cage seal ring (6).
3. Fluid dripping/running out of nozzle when not in use.	Applicator is hanging at end of feed tube when not in use.	 Hang the applicator at same height or higher than the off take point of the feed tube on the container of product. This ensures the delivery valve is free of load which can cause the product to leak past the valve assembly.
4. Hard delivery stroke pressure.	Foreign matter in delivery valve spring or blockage in nozzle.	Remove the nozzle (1 or 2) and ensure the delivery holes are clean and clear. Clean delivery valve (4) and spring (3) and ensure it is positioned in the correct direction.
5. Air is being drawn into the cylinder from in and around the piston.	Foreign matter lodged on or around piston o-ring.	Replace piston o-ring (10 & 12) and lubricate liberally with NJ Phillips Lubricant. <b>DO NOT</b> use petroleum based lubricants – only use synthetic or vegetable oil.
	Piston inlet valve not seated correctly.	Remove the cylinder (7). Hold the rear of the push rod (17) and tighten the piston (11) by firmly holding with a dry cloth and turning clockwise.
	Feed tube not connected to applicator or draw-off correctly.	Ensure the feed tube spring is located correctly on the applicator and draw-off. 
	Feed tube perforated or damaged.	Replace the feed tube.
	Feed tube connection at container or applicator is split or damaged.	Replace container fitting to ensure an air tight seal. Cut feed tube for clean ends.
6. Piston not returning fully on filling stroke.	The piston o-ring and lubricating washer are dry.	Remove the cylinder (7), soak the piston o-ring (10) and lubricating washer (12) in NJ Phillips Lubricant, or synthetic or vegetable oil.
	Blockage in inlet line.	Check the inlet valve (22), feed tube and container draw-off fitting for foreign matter.
	Kinking or restriction of feed tube.	Remove restriction or reposition feed tube. Ensure the springs are in correct position.
	Chemical container not collapsing as applicator draws fluid.	Vent pack or use an NJ Phillips Vented Draw-Off system.

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