

30mL

POUR-ON



njphillips™
DATAMARS




- Administers a variety of animal health pour-on products.
- Easy access for cleaning with external valves and removable cylinder.
- Durable and reliable, with an ergonomic 'ripple-grip' handle.

30mL Pour-On Parts Breakdown

NO SPARE PARTS KITS AVAILABLE.

- | | |
|--------------------------------|-------------------------------|
| 1. Nozzle Tip | 14. Dose Adjuster |
| 2. Valve Spring | 15. 30mL Dose Adjuster Sleeve |
| 3. Valve | 16. Dose Collar |
| 4. Delivery Cage Lock Nut | 17. Dose Adjuster Guide Ring |
| 5. 30mL Barrel F/D | 18. 65mL Ported Pushrod |
| 6. O-Ring (111-S70) | 19. Handle |
| 7. O-Ring (113-S70) | 20. Return Spring |
| 8. Delivery Cage | 21. Lever |
| 9. Washer Piston | 22. Lever Pads |
| 10. 65mL Valve Retainer Piston | 23. Lever Pin (Long) |
| 11. O-Ring (121-F70) | 24. Nozzle Tip (Fan) |
| 12. Felt Washer | 25. Nozzle Tip (Jet) |
| 13. 30-40mL Piston | 26. Feed Tube Springs |



Product	Dose Range	Increments	Material	Packaged Dimensions (WxHxD)	Packaging Type	Code	SAP Code	EAN Barcode	UPC Barcode
30mL Pour-On	2.5mL - 30mL	2.5, 3, 4, 5, 6, 7.5, 8, 9, 10, 11, 12, 12.5, 13, 14, 15, 16, 17, 17.5, 18, 19, 20, 22.5, 25, 27.5, 30mL	Plastic	230 x 380 x 50mm	Headercard & bag	PAM1038	850 0000-627	9 323837 009109	8 12465 01288 3

For more information, email: ahdsinfo@datamars.com or phone toll free (Australia): 1800 247 175

www.njphillips.com